

ENVIKIT BORAX

ENVIKIT 21Bor




- Thank you for choosing products of Binh Lan ENVIBIOCHEM. Please read carefully before proceeding to ensure accurate results.
- Reagents contain acid. In case of contact, immediately flush skin with plenty of water.

Contain: 20 ml of Reagent 1, 50 turmeric strips, a plastic cup, a plastic spoon, a paper-clip and instruction manual.

Application: detection of borax in food such as meat and meat products (minced meat, minced fish, fish cake, meatball, pork, beef, sausage, etc.), cakes, preserved fruits, pickled vegetable, etc. **Sensitivity of Test Kit:** Detection limit: 80 mg/L (80 ppm) **Number of Test / Kit:** 50 Tests

Instruction manual



for 15 minutes. If the strip turns orange or red (a, b), such sample contains borax.

Impact of Borax on Human Health

Consumption of borax-contaminated food is health hazard.

The acute effects are nausea, vomiting and diarrhea, with delayed effects of skin redness and peeling. 15 grams of borax for adults and 5 grams for children can cause hematemesis (blood vomiting) and death.

Accumulation of borax is suspected of damaging fertility or the unborn child. Birth defects include brain and skeletal abnormalities and infertility of male offspring. Animal ingestion studies in several species, at high doses, indicate that borax causes reproductive and developmental effects. According to the U.S. Environmental Protection Agency, animal studies have shown that boric acid has an affinity for male testicles. Poisoning produces inflammation and atrophy of tubules and lowers sperm count. Repeated low dose exposure borax used as a food preservative in noodles, steamed rice, caviar and meatballs, has been implicated in cases of liver cancer over five to ten years.

Borax has been banned as a preservative for food in the United States and many other countries, but they are still used illegally.

Borax Envikit of Binh Lan Envibiochem is applied for detection of Borax in food. The test kit can be easily used outside the laboratory with fast and reliable result.

SAFE YOUR SALVES

BINH LAN ENVIBIOCHEM COMPANY Limited

Web: http://www.binhlan.com Email: binhlanbiochem.gmail.com
Bank account: ACB Bank, Binh Lan Company, STK:
Office 1: Science-Technology Business Incubator of HoChiMinh City 268 Ly Thuong Kiet Street, District 10, HoChiMinh City, Vietnam
Office 2: 595/208A Cach mang thang 8, Ward 15, District 10, HoChiMinh, Vietnam
Tel: 08-39700106, Hotline: 0985714580

For technical advice, price and ordering

Email to binhlanbiochem.gmail.com